



# Euroguard

## Kirurgisk munnbind

Ear loops, blå.

**Art. nr: 17215**

**Forpakning**

50/boks, 500 per kartong

Indre boks: 842942003868

Ytre boks: 10842942003865

**CE-merket**

Spesifikasjoner	
EN14683 klassifisering	Type II
BFE	> 98% (at 3.0 microns)
PFE	> 98% (at 0.1 microns)
Delta P	22.5 Pa
Væskemotstand	80 mm Hg
Brennbarhet	Meets Class 1 (16 CFR 1610)
Dimensjoner	17.5 x 9.3 cm
Ytre lag	Polyspun bond polypropylene
Filter lag	Melt blown polypropylene
Innvendig fôr	Polyspun bond polypropylene
Rør (piping)	Spunlace (viscose/PET) top. Self piping on bottom
Konstruksjon	All seams ultrasonically welded – Total 3 layers
Nese-stykke	Aluminum alloy strip, 14.0 x 0.3 cm
Ear-loops	Spandex / polyester Rounded type 13.0 cm free length Ultrasonically welded to outer side

Testing: Performed at Nelson Labs, Salt Lake City (USA) in accordance with ASTM F2100-07 and EN14683 standards